

Accent™ SR

Single-Chamber Pacemaker

MODEL PM1110



SPECIFICATIONS

- A two-tone audible alert allows programming for the patient to be alerted to changes in device performance, or information can be sent directly to the clinician through the Merlin.net™ Patient Care Network (PCN).
- AutoCapture™ Pacing System offers the maximum in threshold adaptability and patient safety with ventricular Beat-by-Beat™ capture confirmation. The AutoCapture Pacing System automatically delivers a 5,0 V backup safety pulse when noncapture is detected, and it may be programmed to either a bipolar or unipolar configuration.
- Exclusive SenseAbility™ feature, with Decay Delay and Threshold Start, provides the flexibility to fine-tune sensing to individual patient needs and help eliminate oversensing of T waves, fractionated QRS complexes, and other extraneous signals.
- Real-time EGM waveform, as well as the associated event markers that precede and follow a specific triggering event, can be programmed to automatically record up to 14 minutes of stored EGMs when encountering one or more programmable trigger options.
 - Preferential EGM storage capability allows prioritisation of episode storage.
- Weekly Lead Impedance Trend displays the current measurement, historical test results, pacing polarity and any polarity switches.
- Upon interrogation, the device displays the last automatically measured capture threshold results from the ventricle. In addition, the pacemaker automatically measures intrinsic R-wave activity daily and displays the last test results in combination with a weekly R-wave trend. Results are displayed with follow-up EGMs for quick verification.

Indications: Implantation is indicated in one or more of the following permanent conditions: syncope, presyncope, fatigue, disorientation due to arrhythmia/bradycardia, or any combination of those symptoms. **Rate-Modulated Pacing** is indicated for patients with chronotropic incompetence, and for those who would benefit from increased stimulation rates concurrent with physical activity. **Atrial Pacing** is indicated for patients with sinus node dysfunction and normal AV and intraventricular conduction systems. **Ventricular Pacing** is indicated for patients with significant bradycardia and normal sinus rhythm with only rare episodes of A-V block or sinus arrest, chronic atrial fibrillation, severe physical disability.

Contraindications:

Implanted Cardioverter-Defibrillator (ICD). Devices are contraindicated in patients with an implanted cardioverter-defibrillator. **Rate-Adaptive Pacing** may be inappropriate for patients who experience angina or other symptoms of myocardial dysfunction at higher sensor-driven rates. An appropriate Maximum Sensor Rate should be selected based on assessment of the highest stimulation rate tolerated by the patient. **Single-Chamber Ventricular Demand Pacing** is relatively contraindicated in patients who have demonstrated pacemaker syndrome, have retrograde VA conduction, or suffer a drop in arterial blood pressure with the onset of ventricular pacing. **Single-Chamber Atrial Pacing** is relatively contraindicated in patients who have demonstrated compromise of AV conduction. For specific contraindications associated with individual modes, refer to the programmer's on-screen help.

Warnings and Precautions: To prevent permanent damage to the device and tissue damage at the electrode/tissue interface:

- Electrosurgery. Do not use electrosurgical devices in the vicinity of an implanted device. If electrocautery is necessary, use a bipolar cauteriser or place the indifferent electrode as far from the device as possible.
- Lithotripsy. Do not focus a lithotripsy beam within 16 cm of the device. Program the device to Sensor Off prior to lithotripsy to prevent inappropriate increases in pacing rate. A thorough assessment of device function with special attention to the sensor should be performed following exposure to lithotripsy.
- Therapeutic Radiation. Do not use ionising radiation in the vicinity of an implanted device. Radiation therapy may damage the electronic circuitry of the device.
- Ultrasound Treatment. Do not use therapeutic ultrasound within 16 cm of the device.

Perform a thorough assessment of device function following exposure to any of the above.

Device Communication. Communication with the device can be affected by electrical interference and strong magnetic fields. If this is a problem, turn off nearby electrical equipment or move it away from the patient and the programmer. If the problem persists, contact St. Jude Medical.

External Defibrillation. The electronic circuitry in the device provides protection from defibrillation discharges. Nevertheless, do not place defibrillator paddles directly over the device or pacing lead. Following defibrillation, ensure that the device is operating correctly.

Magnetic Resonance Imaging (MRI). MRI for patients with implantable devices has been contraindicated by MRI manufacturers. Clinicians should carefully weigh the decisions to use MRI with pacemaker patients. Additional safety concerns include:

- Magnetic and RF fields produced by MRI may increase pacing rate, inhibit pacing, cause asynchronous pacing or result in pacing at random rates
- MRI may result in changes in capture thresholds due to heating of pacing leads
- MRI may irreversibly damage the device
- Patients should be closely monitored during the MRI
- Assess the device function before and after exposure to MRI.

CT Scans. CT scans, due to their increased power levels and long exposure times, have the remote possibility of interfering with implanted devices. The potential interference is transient and occurs only when the X-ray signal is present. Continuous exposure may cause a temporary sensor rate increase. In addition, there is a remote possibility for a device to intermittently oversense while the CT scanning beam is directly over the implanted device.

Potential Adverse Events: The following are potential complications associated with the use of any pacing system: arrhythmia, heart block, thrombosis, threshold elevation, valve damage, pneumothorax, myopotential sensing, vessel damage, air embolism, body rejection phenomena, cardiac tamponade or perforation, formation of fibrotic tissue/local tissue reaction, inability to interrogate or program a device because of programmer malfunction, infection, interruption of desired device function due to electrical interference, loss of desired pacing and/or sensing due to lead displacement, body reaction at electrode interface, or lead malfunction (fracture or damage to insulation), loss of normal device function due to battery failure or component malfunction, device migration, pocket erosion, or hematoma, pectoral muscle stimulation, phrenic nerve or diaphragmatic stimulation. The following, in addition to the above, are potential complications associated with the use of rate-modulated pacing systems: inappropriate, rapid pacing rates due to sensor failure or to the detection of signals other than patient activity, loss of activity-response due to sensor failure, palpitations with high-rate pacing.

Refer to the User's Manual for detailed indications, contraindications, warnings, precautions and potential adverse events.



ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.

PHYSICAL SPECIFICATIONS

Models	PM1110
Telemetry	Inductive
Dimensions (mm)	42 x 52 x 6
Weight (g)	18
Volume (cc)	9,5 ¹
Connector	IS-1

PARAMETER SETTINGS

Rate/Timing

Ventricular Pace/Sense Refractory (Fixed) (ms)	125; 160-400 in steps of 30; 440; 470 ²
Base Rate (min ⁻¹)	30-130 in steps of 5; 140-170 in steps of 10
Mode	V00(R); VVI(R); VVT(R); Pacing Off
Hysteresis Rate (min ⁻¹)	Off; 30 ¹ -150 in steps of 5
<i>Search Interval (min⁻¹)</i>	Off; 1; 5; 10; 15; 30
<i>Cycle Count</i>	1-16 by 1
<i>Intervention Rate (min⁻¹)</i>	Off; 80-120 in steps of 10; Intrinsic +0; Intrinsic +10; Intrinsic +20; Intrinsic +30; Same as Base Rate
<i>Intervention Duration (min)</i>	1-10 in 1 minute intervals
<i>Recovery Time</i>	Fast; Medium; Slow; Very Slow
Rest Rate (min ⁻¹)	Off; 30-150; in steps of 5

Output/Sensing

V Pulse Amplitude (V)	0,25-4.0 in steps of 0,25; 4,5-7,5 in steps of 0,5
V Pulse Width (ms)	0,05; 0,1-1,5 in steps of 0,1
V Sensitivity (mV)	0,5-5,0 in steps of 0,5; 6-10 in steps of 1,0; 12,5 ⁴
V Pulse Configuration	Unipolar (tip-case); Bipolar (tip-ring)
V Sense Configuration	Unipolar Tip (tip-case); Bipolar (tip-ring); Unipolar Ring (ring-case)
Ventricular AutoCapture™	
Pacing System	On; Off
<i>Primary Pulse Configuration</i>	Unipolar; Bipolar
<i>Backup Pulse Configuration</i>	Unipolar; Bipolar
<i>Backup Pulse Amplitude (V)</i>	5,0 ⁵
<i>Search Interval (hours)</i>	8; 24
SenseAbility™ Technology	
<i>Max Sensitivity (mV)</i>	Off; On (Automatic Sensitivity Control adjustment for ventricular events)
<i>Threshold Start</i>	(Ventricular Post-Sense) 50; 62,5; 75; 100% (Ventricular Post-Pace) Auto; 0,2-3,0 in steps of 0,1 mV
<i>Decay Delay (ms)</i>	(Ventricular Post-Sense) 0; 30; 60; 95; 125; 160; 190; 220 (Ventricular Post-Pace) Auto; 0; 30; 60; 95; 125; 160; 190; 220

Rate-Modulated Parameters

Maximum Sensor Rate (min ⁻¹)	80-150 in steps of 5; 160-180 in steps of 10
Rate Responsive VREF	Off; Low; Medium; High
Shortest VREF	125-475 in steps of 25
Reaction Time	Very Fast; Fast; Medium; Slow
Recovery Time	Fast; Medium; Slow; Very Slow
Sensor	On; Off; Passive
Slope	Auto (-1); Auto (+0); Auto (+1); Auto (+2); Auto (+3); 1-16 in steps of 1
Threshold	Auto (-0,5); Auto (+0,0); Auto (+0,5); Auto (+1,0); Auto (+1,5); Auto (+2,0); 1-7 in steps of 0,5

Stored Electrograms

Options	
<i>Priority Options</i>	Off; Low; High
<i>Channel</i>	1; 2; 3
Triggers	
<i>Magnet Response</i>	Off; Low; High
<i>High Ventricular Rate Rate (min⁻¹)</i>	Off; Low; High
<i>No. of Consecutive Cycles</i>	125-300 in steps of 25
<i>Advanced Hysteresis</i>	2; 3; 4; 5; 10; 15; 20
<i>Noise Reversion</i>	Off; Low; High

Other

Lead Monitoring	Monitor; Auto Polarity Switch
<i>V Low Impedance Limit (Ω)</i>	100-500 in steps of 25
<i>V High Impedance Limit (Ω)</i>	750-2500 in steps of 250; 3000
Magnet Response	Off; Battery Test
Lead Type	Uncoded; Unipolar; Bipolar
NIPS Options	
<i>Stimulation Chamber</i>	Ventricular
<i>Coupling Interval (ms)</i>	100-800 in steps of 10
<i>S1 Count</i>	2-25 in steps of 1
<i>S1²; S2; S3 and S4 Cycle (ms)</i>	Off; 100-800 in steps of 10 (Fixed or Adaptive)
Diagnostic Trends	Exercise; Lead Impedance; R Wave; V Threshold

Patient Notifiers

Programmable Notifiers (On; Off)	Device at ERI; Ventricular Lead Impedance Out of Range
Device Reset	On
Entry into Backup VVI Mode	On
Audible Duration (sec)	2; 4; 6; 8; 10; 12; 14; 16
Number of Audible Alerts per Notification	2
Number of Notifications	1-16
Time Between Notifications (hours)	10; 22

- 1 ± 0,5 cc
- 2 Programming options dependent on pacing mode.
- 3 The highest available setting for hysteresis rate will be 5 min⁻¹ below the programmed base rate.
- 4 Sensitivity is with respect to a 20 ms haversine test signal.
- 5 This parameter is not programmable.
- 6 S1 Burst Cycle is applied at the preprogrammed S1 cycle length.

ATRIAL FIBRILLATION CARDIAC RHYTHM MANAGEMENT CARDIOVASCULAR NEUROMODULATION

Global Headquarters
One Lillehei Plaza
St. Paul, Minnesota 55117
USA
+1 651 483 2000
+1 651 490 4310 Fax

Cardiac Rhythm Management Division
15900 Valley View Court
Sylmar, California 91342
USA
+1 818 362 6822
+1 818 364 5814 Fax

St. Jude Medical Sweden AB
Veddestavägen 19
175 84 Järfälla
Sweden
+46 8 474 40 00
+46 8 760 95 42 Fax

St. Jude Medical Europe, Inc.
The Corporate Village
Figueras Building
Avenue Da Vinci laan, 11 Box F1
B-1935 Zaventem
Belgium
+32 2 774 68 11
+32 2 772 83 84 Fax

St. Jude Medical Brasil Ltda.
Rua Frei Caneca, 1380
7º ao 9º andares
01307-002 - São Paulo (SP)
Brazil
+55 11 5080 5400
+55 11 5080 5423 Fax

St. Jude Medical (Hong Kong) Ltd.
Unit 2701-07 27/F, COSCO Tower
Grand Millennium Plaza
183 Queen's Road
Central, Hong Kong
+852 2996 7688
+852 2956 0622 Fax

St. Jude Medical Japan Co., Ltd.
3-1-30, Minami-Aoyama
Minato-ku
Tokyo 107 0062
Japan
+81 3 3423 6450
+81 3 3402 5586 Fax

SJMprofessional.com



Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country. Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK. are trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2009 St. Jude Medical, Inc. All Rights Reserved.
Item GMC454